

SGK2 Polyclonal Antibody

Catalog number: 11185-1-AP

Size: 45 µg/150 µl


Source: Rabbit

Isotype: IgG


Synonyms:

SGK2; dJ138B7.2, H SGK2,

SGK2



Immunohistochemical of paraffin-embedded human pancreas cancer using 11185-1-AP (SGK2 antibody) at dilution of 1:50 (under 10x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with 11185-1-AP (SGK2 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, IHC, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated SGK2 MW:	48 kDa
Observed SGK2 MW:	
Positive WB detected in	Transfected HEK-293 cells
Positive IHC detected in	Human pancreas cancer tissue
Recommended dilution:	WB: 1:500-1:5000
	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1669
GenBank accession number:	BC014037
Gene ID (NCBI):	10110
Full name:	Serum/glucocorticoid regulated kinase 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.