

NUP93 Polyclonal Antibody

Catalog number: 13077-1-AP


Size: 33 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

NUP93; KIAA0095; MGC21106



HEK-293 cells were subjected to SDS PAGE followed by western blot with 13077-1-AP(NUP93 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HEK-293 cells, A431 cells
Calculated NUP93 MW:	819aa,93kd
Observed NUP93 MW:	93kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3735
GenBank accession number:	BC034346
Gene ID (NCBI):	9688
Full name:	Nucleoporin 93kDa

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.