

## MUC20 Polyclonal Antibody

Catalog number: 13380-1-AP

Size: 47 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

MUC20; FLJ14408; FLJ53153;

KIAA1359



human liver tissue were subjected to SDS PAGE followed by western blot with 13380-1-AP(MUC20 Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Positive WB detected in:	Human liver tissue, human kidney tissue
Calculated MUC20 MW:	503aa,55kd
Observed MUC20 MW:	70-75kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag4205
GenBank accession number:	BC029267
Gene ID (NCBI):	200958
Full name:	Mucin 20, cell surface associated

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.