

**SAFETY DATA SHEET**

according to Regulation (EC) No. 1907/2006

**AccuGENE® 10X TAE Buffer**

Version 2.0

Revision Date 07.12.2018

Print Date 14.01.2020

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1 Product identifier**

Product name : AccuGENE® 10X TAE Buffer  
Material number : 50841  
Synonyms : Accugene 10X TAE Buffer

**1.2 Relevant identified uses of the substance or mixture and uses advised against**

Use of the Substance/Mixture : Scientific research and development

**1.3 Details of the supplier of the safety data sheet**

Company : **Lonza Ltd**  
Muenchensteinerstrasse 38  
CH-4002 Basel, Switzerland  
Business Telephone: +41 61 316 81 11

**Lonza Verviers Sprl**  
Parc Industriel de Petit-Rechain  
BE-4800 Verviers, Belgium  
Business Telephone: +32 8732 1611

**Lonza Cologne GmbH**  
Nattermannallee 1  
DE-50829 Köln, Germany  
Business Telephone: + 49 221 99 1990

**Lonza Copenhagen ApS**  
Strandhaven 12  
DK-2665 Vallensbaek Strand, Denmark  
Business Telephone: + 45 4356 7400

E-mail address : sds@lonza.com

Responsible/issuing person

**1.4 Emergency telephone number**

Emergency telephone number : Lonza Ltd, CH-4002 Basel, Switzerland  
Telephone: +41 61 313 94 94 (24h)

**SECTION 2: Hazards identification****2.1 Classification of the substance or mixture****Classification (REGULATION (EC) No 1272/2008)**

Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.

**2.2 Label elements****Labelling (REGULATION (EC) No 1272/2008)**

Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.

## AccuGENE® 10X TAE Buffer

### Additional Labelling

EUH210 Safety data sheet available on request.

### 2.3 Other hazards

This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

## SECTION 3: Composition/information on ingredients

### 3.2 Mixtures

#### Components

Chemical name	CAS-No. EC-No. Index-No. REACH Registration Number	Classification	Concentration (% w/w)
Acetic acid	64-19-7  607-002-00-6 01-2119475328-30-XX XX	Flam. Liq. 3; H226 Met. Corr. 1; H290 Acute Tox. 4; H312 Skin Corr. 1A; H314 Eye Dam. 1; H318	>= 1 - < 3

For explanation of abbreviations see section 16.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

- If inhaled : Remove to fresh air.  
If breathing is irregular or stopped, administer artificial respiration.
- In case of skin contact : After contact with skin, wash immediately with plenty of soap and water.  
Get medical attention if irritation develops and persists.
- In case of eye contact : In the case of contact with eyes, rinse immediately with plenty of water and seek medical advice.  
Call a physician if irritation develops or persists.
- If swallowed : Call a physician immediately.  
Clean mouth with water and drink afterwards plenty of water.  
Do not induce vomiting without medical advice.  
Never give anything by mouth to an unconscious person.

### 4.2 Most important symptoms and effects, both acute and delayed

Symptoms : No information available.

### 4.3 Indication of any immediate medical attention and special treatment needed

Treatment : No information available.

**AccuGENE® 10X TAE Buffer**

---

**SECTION 5: Firefighting measures****5.1 Extinguishing media**

Suitable extinguishing media : Dry powder  
Water spray  
Foam

**5.2 Special hazards arising from the substance or mixture**

Specific hazards during firefighting : Heating or fire can release toxic gas.

**5.3 Advice for firefighters**

Special protective equipment for firefighters : In the event of fire, wear self-contained breathing apparatus.

Further information : Use water spray to cool unopened containers.

---

**SECTION 6: Accidental release measures****6.1 Personal precautions, protective equipment and emergency procedures**

Personal precautions : Use respirator when performing operations involving potential exposure to vapour of the product.

**6.2 Environmental precautions**

Environmental precautions : Prevent product from entering drains.

**6.3 Methods and material for containment and cleaning up**

Methods for cleaning up : Do not flush into surface water or sanitary sewer system.  
Soak up with inert absorbent material.

**6.4 Reference to other sections**

For personal protection see section 8.

---

**SECTION 7: Handling and storage****7.1 Precautions for safe handling**

Advice on safe handling : Provide sufficient air exchange and/or exhaust in work rooms.  
Avoid contact with skin and eyes.

Advice on protection against fire and explosion : Take precautionary measures against static discharges.

Hygiene measures : Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product.

**7.2 Conditions for safe storage, including any incompatibilities**

Requirements for storage areas and containers : Keep container closed when not in use.

**7.3 Specific end use(s)**

## AccuGENE® 10X TAE Buffer

Specific use(s) : No information available.

### SECTION 8: Exposure controls/personal protection

#### 8.1 Control parameters

Countries not listed may have their own country specific values.

#### Occupational Exposure Limits

##### European Union

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	STEL	20 ppm 50 mg/m <sup>3</sup>	ECTLV

##### Switzerland

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	TWA	10 ppm 25 mg/m <sup>3</sup>	SMAK
		STEL	20 ppm 50 mg/m <sup>3</sup>	SMAK

##### Germany

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	MAK	10 ppm 25 mg/m <sup>3</sup>	DFG MAK
Peak-limit: excursion factor (category)		Peak-limit: excursion factor (category) 2		
			10 ppm 25 mg/m <sup>3</sup>	TRGS 900
Peak-limit: excursion factor (category)		Peak-limit: excursion factor (category) 2		

##### Belgium

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	STEL	15 ppm 38 mg/m <sup>3</sup>	BE/OEL
		TWA	10 ppm 25 mg/m <sup>3</sup>	BE/OEL

##### Denmark

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	GV	10 ppm 25 mg/m <sup>3</sup>	GV (DK)

##### France

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	VLE	10 ppm 25 mg/m <sup>3</sup>	FVL

##### Italy

## AccuGENE® 10X TAE Buffer

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	TWA	10 ppm 25 mg/m <sup>3</sup>	OEL (IT)
		STEL	20 ppm 50 mg/m <sup>3</sup>	OEL (IT)

### Netherlands

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	TGG	25 mg/m <sup>3</sup>	NL OEL

### Spain

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	STEL	15 ppm 37 mg/m <sup>3</sup>	VLA
		TWA	10 ppm 25 mg/m <sup>3</sup>	VLA

### Sweden

Components	CAS-No.	Value type (Form of exposure)	Control parameters	Basis
Acetic acid	64-19-7	TGV	10 ppm 25 mg/m <sup>3</sup>	SWO
		TWA	5 ppm 13 mg/m <sup>3</sup>	SWO

### United Kingdom

#### Derived No Effect Level (DNEL) according to Regulation (EC) No. 1907/2006:

Substance name	End Use	Exposure routes	Potential health effects	Value	
Ethylenediaminetetraacetic acid disodium salt dihydrate	Workers	Inhalation	Acute local effects	3 mg/m <sup>3</sup>	
	Consumers	Oral	Long-term systemic effects	25 mg/kg	
	Consumers	Inhalation	Long-term local effects	0,6 mg/m <sup>3</sup>	
	Consumers	Inhalation	Acute local effects	1,2 mg/m <sup>3</sup>	
	Workers	Inhalation	Long-term local effects	1,5 mg/m <sup>3</sup>	
	Acetic acid	Workers	Inhalation	Short-term exposure, Local effects	25 mg/m <sup>3</sup>
	Workers	Inhalation	Long-term local effects	25 mg/m <sup>3</sup>	
	Consumers	Inhalation	Long-term local effects	25 mg/m <sup>3</sup>	
	Consumers	Inhalation	Short-term exposure, Local effects	25 mg/m <sup>3</sup>	
	Workers	Skin contact	Long-term local effects	10 mg/m <sup>3</sup>	

#### Predicted No Effect Concentration (PNEC) according to Regulation (EC) No. 1907/2006:

Substance name	Environmental Compartment	Value
Ethylenediaminetetraacetic acid disodium salt dihydrate	Fresh water	2,2 mg/l
	Soil	0,72 mg/kg
	Sewage treatment plant	43 mg/l
	Intermittent use/release	1,2 mg/l
	Marine water	0,22 mg/l
	Acetic acid	Fresh water sediment
	Sewage treatment plant	85 mg/l

## AccuGENE® 10X TAE Buffer

	Soil	0,478 mg/kg
	Intermittent use/release	30,58 mg/l
	Fresh water	3,058 mg/l
	Marine water	0,3058 mg/l
	Marine sediment	1,136 mg/kg

### 8.2 Exposure controls

#### Engineering measures

Use only in area provided with appropriate exhaust ventilation.

#### Personal protective equipment

Eye protection : Tightly fitting safety goggles

Hand protection

Material : Nitrile rubber

Rate of permeability : > 480 min

Skin and body protection

: Choose body protection according to the amount and concentration of the dangerous substance at the work place. No special protective equipment required.

Respiratory protection

: In the case of vapour formation use a respirator with an approved filter.

Respirator with a vapour filter (EN 141)

Respirator with ABEK filter.

## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

Appearance : Aqueous solution

Colour : colourless, clear

Odour : odourless

Odour Threshold : no data available

pH : 7 - 9

Melting point/freezing point : no data available

Boiling point/boiling range : ca. 100 °C

Flash point : no data available  
Water

Evaporation rate : no data available

Flammability (solid, gas) : no data available

Upper explosion limit : no data available

Lower explosion limit : no data available

Vapour pressure : no data available

## AccuGENE® 10X TAE Buffer

---

Relative vapour density	:	no data available
Relative density	:	no data available
Density	:	1,02 g/cm <sup>3</sup>
Solubility(ies)		
Water solubility	:	soluble
Partition coefficient: n-octanol/water	:	no data available
Auto-ignition temperature	:	no data available
Decomposition temperature	:	no data available
Viscosity	:	no data available
Explosive properties	:	no data available
Oxidizing properties	:	no data available

### 9.2 Other information

no data available

---

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

Stable under recommended storage conditions.

### 10.2 Chemical stability

Stable under normal conditions.

### 10.3 Possibility of hazardous reactions

Hazardous reactions : Stable under normal conditions.

### 10.4 Conditions to avoid

Conditions to avoid : no data available

### 10.5 Incompatible materials

Materials to avoid : Strong oxidizing agents

### 10.6 Hazardous decomposition products

No decomposition if stored normally.

---

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

Information on likely routes of exposure : Eyes  
Ingestion

#### Acute toxicity

Acute dermal toxicity : Acute toxicity estimate: > 2 000 mg/kg  
Method: Calculation method

## AccuGENE® 10X TAE Buffer

---

### Further information

Remarks: No data is available on the product itself.

Remarks: Health injuries are not known or expected under normal use.

---

## SECTION 12: Ecological information

### 12.1 Toxicity

Toxicity to fish : Remarks: no data available

### 12.2 Persistence and degradability

Biodegradability : Result: no data available

### 12.3 Bioaccumulative potential

Bioaccumulation : Remarks: no data available

### 12.4 Mobility in soil

Distribution among environmental compartments : Remarks: no data available

### 12.5 Results of PBT and vPvB assessment

Assessment : This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

### 12.6 Other adverse effects

Additional ecological information : no data available

---

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

Product : Contact waste disposal services.  
Dispose of in accordance with local regulations.

Contaminated packaging : Dispose of as unused product.

**AccuGENE® 10X TAE Buffer**

---

**SECTION 14: Transport information**

**IATA** Not dangerous goods

- 14.1 UN number : Not applicable
- 14.2 Proper shipping name : Not applicable
- 14.3 Transport hazard class : Not applicable
- 14.4 Packing group : Not applicable
- 14.5 Environmental hazards : no

**IMDG** Not dangerous goods

- 14.1 UN number : Not applicable
- 14.2 Proper shipping name : Not applicable
- 14.3 Transport hazard class : Not applicable
- 14.4 Packing group : Not applicable
- 14.5 Environmental hazards : Marine pollutant: no

**ADR** Not dangerous goods

- 14.1 UN number : Not applicable
- 14.2 Proper shipping name : Not applicable
- 14.3 Transport hazard class : Not applicable
- 14.4 Packing group : Not applicable
- 14.5 Environmental hazards : no

**RID** Not dangerous goods

- 14.1 UN number : Not applicable
- 14.2 Proper shipping name : Not applicable
- 14.3 Transport hazard class : Not applicable
- 14.4 Packing group : Not applicable
- 14.5 Environmental hazards : no

**DOT** Not dangerous goods

- 14.1 UN number : Not applicable
- 14.2 Proper shipping name : Not applicable
- 14.3 Transport hazard class : Not applicable
- 14.4 Packing group : Not applicable
- 14.5 Environmental hazards : no

## AccuGENE® 10X TAE Buffer

---

**TDG** : Not dangerous goods

**14.1 UN number** : Not applicable

**14.2 Proper shipping name** : Not applicable

**14.3 Transport hazard class** : Not applicable

**14.4 Packing group** : Not applicable

**14.5 Environmental hazards** : no

**14.6 Special precautions for user** : none

**14.7 Transport in bulk according to Annex II of Marpol and the IBC Code** : Not applicable

---

### SECTION 15: Regulatory information

#### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

REACH - Restrictions on the manufacture, placing on the market and use of certain dangerous substances, preparations and articles (Annex XVII) : Not applicable

REACH - Candidate List of Substances of Very High Concern for Authorisation (Article 59). : Not applicable

REACH - List of substances subject to authorisation (Annex XIV) : Not applicable

Regulation (EC) No 1005/2009 on substances that deplete the ozone layer : Not applicable

Regulation (EC) No 850/2004 on persistent organic pollutants : Not applicable

Seveso III: Directive 2012/18/EU of the European Parliament and of the Council on the control of major-accident hazards involving dangerous substances.  
Not applicable

#### National regulatory information

Water contaminating class (Germany) : WGK 1 slightly hazardous to water  
Classification according to AwSV, Annex 1 (5.2)

#### 15.2 Chemical safety assessment

not required

---

### SECTION 16: Other information

#### Full text of H-Statements

H226 : Flammable liquid and vapour.

H290 : May be corrosive to metals.

H312 : Harmful in contact with skin.

H314 : Causes severe skin burns and eye damage.

H318 : Causes serious eye damage.

## AccuGENE® 10X TAE Buffer

---

### Further information

#### Full text of other abbreviations

ADR - European Agreement concerning the International Carriage of Dangerous Goods by Road; ATE - Acute Toxicity Estimate; AwSV - Verordnung über Anlagen zum Umgang mit wassergefährdenden Stoffen / Ordinance on facilities for handling substances that are hazardous to water; BPR – Biocidal Product Regulation; bw - Body weight; CAS - Chemical Abstract Service; CLP - Classification Labelling Packaging Regulation, Regulation (EC) No 1272/2008; CMR - Carcinogen, Mutagen or Reproductive Toxicant; DIN - Standard of the German Institute for Standardisation; DNEL-Derived No Effect Level; DOT - Department of Transportation; EC – European Community; ECHA - European Chemicals Agency; EC-Number - European Community number; ECx - Concentration associated with x% response; EmS (Emergency Response Procedures for Ships Carrying Dangerous Goods); EN – European Standard; ErCx - Concentration associated with x% growth rate response; ERG - Emergency Response Guide; EU OEL - European Occupational Exposure Limit; GHS -Globally Harmonized System of Classification and Labelling of Chemicals; GLP - Good Laboratory Practice; GV – Danish Exposure Limits for Substances and Materials; IARC - International Agency for Research on Cancer; IATA - International Air Transport Association; IBC - International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk; IC50 - Half maximal inhibitory concentration; IMDG - International Maritime Dangerous Goods; LC50 - Lethal Concentration to 50 % of a test population; LD50 - Lethal Dose to 50% of a test population (Median Lethal Dose); LOFT - Danish Threshold Limit Value; MAK - German Threshold Limit Value; MARPOL - International Convention for the Prevention of Pollution from Ships; n.o.s. - Not Otherwise Specified; NIOSH/Guide – National Institute of Safety and Health Guidebook; NO(A)EC - No Observed (Adverse) Effect Concentration; NO(A)EL - No Observed (Adverse) Effect Level; NTP - National Toxicology Program; OECD - Organization for Economic Co-operation and Development; PBT - Persistent, Bioaccumulative and Toxic substance; PEL - Permissible Exposure Limit; PNEC - Predicted no Effect Concentration; REACH - Regulation (EC) No 1907/2006 of the European Parliament and of the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals; REL - Recommended Exposure Limit; RID - Regulations concerning the International Carriage of Dangerous Goods by Rail; STEL - Short-Term Exposure Limit; TDG - Transportation of Dangerous Goods; TGG – Dutch Threshold Limit Value; TGV – Swedish OEL; TLV Threshold Limit Value; TLV-C - Threshold Limit Value Ceiling; TWA -Time Weighted Average; UDS - Unscheduled DNA Synthesis; UN - United Nations; VLE - Valeurs limites d'exposition professionnelle aux agents chimiques en France; VME - Valeur (Limite) Moyenne d'Exposition; VOC - Volatile Organic Compound[s]; WEEL - Workplace Environmental Exposure Level; % w/w Percent weight by weight; %(V) Percent Volume

Date format : dd.mm.yyyy

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

99 / EN