

空运
By Air



NO.2020034801



货物运输条件鉴定书

Certification
for Safe Transport of Chemical Goods

危险品

样品名称： 一次性采样管

Sample Name: Disposable Sampler

委托单位： 无锡耐思生物科技有限公司

生产单位： 无锡耐思生物科技有限公司



上海化工院检测有限公司

Shanghai Research Institute of Chemical Industry Testing Co., Ltd

空运
By Air

空运
By Air

货物运输条件鉴定书

Certification for Safe Transport of Chemical Goods

NO. 2020034801

Page 1 / 2

样品名称 Sample Name	中文 Chinese	一次性采样管
	英文 English	Disposable Sampler
委托单位 Consignor	无锡耐思生物科技有限公司	
生产单位 Manufacturer	无锡耐思生物科技有限公司	
检验方法、程序 Inspection Methods and Procedures	国际航空运输协会《危险品规则》61版 IATA Dangerous Goods Regulations (DGR) 61st Edition	
样品外观与气味 Appearance & Odor	类白色半透明塑料外壳（内含无色透明液体），稍有气味 Off-white semi-transparent Plastic shell (containing colorless transparent liquid), Weak odor	
I D E N T I F I C A T I O N 鉴 定 结 论 C O N C L U S I O N	1. 危险性识别 (Hazards identification) 易燃液体。 Flammable liquid.	
	2. 空运按照 IATA DGR 办理的类型 (Suggestion according to IATA DGR) Shipping Name: Ethanol solution Class or Division: 3 UN Number: UN1170	
	3. 包装要求 (Packaging requirements) 按 III 类包装要求办理。 Packing Group III.	
	检验日期: Inspection Date: 2020-03-09 签发日期: Issue Date: 2020-03-09 生效日期: Effective Date: 2020-03-09	
备注 Comment	无。 None.	

批准
Approver: 张一明

审核
Checker: 董学胜

主检
Appraiser: 周敏



货物运输条件鉴定书

Certification for Safe Transport of Chemical Goods

NO. 2020034801

Page 2 / 2

鉴定项目 Identification Items	鉴定结果 Identification Conclusion Results
爆炸危险性鉴定 Identification of Explosive Hazard	该货物不属于爆炸品。 The product is not classified in Explosives.
易燃危险性鉴定 Identification of Flammable Hazards	经测试, 闭杯闪点51度, 表明该货物属于3类易燃液体。 In the closed-cup flash point test, fp 51°C, so the substance is classified in Class 3 (Flammable Liquids).
氧化危险性鉴定 Identification of Oxidative Hazards	该货物不属于氧化剂和有机过氧化物。 The product is not classified in oxidizing substances and organic peroxides.
毒害及传染危险性鉴定 Identification of Toxic & Infectious Hazards	该货物不属于有毒和感染性物质。 The product is not classified in toxic and infectious substances.
放射危险性鉴定 Identification of Radioactive Hazard	该货物无放射危险性。 The product is not classified in radioactive material.
腐蚀危险性鉴定 Identification of Corrosive Hazard	该货物不属于腐蚀品。 The product is not classified in corrosives.
其他危险性鉴定 Identification of other Hazards	该货物无其它危险性。 The product presents no other dangerous properties.

-验证码:542992-

报告结束

