

RGS10 Polyclonal Antibody

Catalog number: 11163-1-AP


Size: 25 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

RGS10; RGS10



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 11163-1-AP(RGS10 antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated RGS10 MW:	21 kDa
Observed RGS10 MW:	20-22kd
Positive WB detected in	Mouse thymus tissue, K-562 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1644
GenBank accession number:	BC009361
Gene ID (NCBI):	6001
Full name:	Regulator of G-protein signaling 10

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.