



# Large Sky Line Digital Orbital Shaker DOS-10L / 20L



## Product Description

The Sky Line series is a high-performance line of orbital and rocking shakers capable of non-stop operation in a wide variety of medical, chemical, and biological laboratories. Designed for flexibility, these units come with a great selection of optional accessories and attachments, and in an assortment of sizes and loading capacities to satisfy your specific lab needs. With an excellent record of reliability and long life, the Sky Line series will be a superior partner in your laboratory.

## Features

- Optimal amplitude of platform rotation for intensive mixing
- Advanced stepper motor control system
- Digital timer
- High-precision, rotational speed controls
- Sound notification when mixing complete
- Onboard memory for storing settings

## Specifications

Size of Working Space of Platform	415 x 295
Amplitude of Rotation of the Platform (Rocking Angle)	10 / 20
Type of Rotation	Orbital
Time Controller	from 1 to 599 / limitless
Speed of Rotation (Rocking)	from 1 to 500 / 300
Environmental Temperature	from +4 to +45
Relative Moisture of Atmosphere (at the temperature of 20 °C)	up to 80
Power Supply (Adapter)	DC 24V, 1A
Size of the Device (length x width x height)	415 x 365 x 100
Weight	7.7/8.7

