

Product Description: safeForce fluid resistant isolation gowns are constructed with 35g polyethylene coated polypropylene material. Ultrasonic welded seams provide added protection from liquid splash. Elastic cuffs help ensure a secure fit and extra long ties are sewn at the neckline and waist for convenience. Open back allows air circulation for a cooler, more comfortable fit.


Product Size: XL, 47" Length (120 x 140cm)

Packaging: 10 per Bag; 10 Bags/CS; 100/CS

Colors: Yellow (20A00M242), Blue (20A00M243)

Material Description: PP+PE laminated 35g

- Features:**
- Fluid Resistant
 - Ultrasonically Welded Seams
 - Low Linting
 - Complies with ANSI/AAMI PB70I Level 2 Requirements



Measurements	Dimensions (Unit: cm)									
	A	B	C	E	F	H	J	M	P	G
XL 120x140	120	140	60	24	27	62	6	65	110	60
Tolerance	±2	±2	±1	±1	±1	±2	±0.5	±2	±2	±2