

Product Name:	ABS SSC 34
Description:	<p> ABS SSC Class Includes 6 Canisters Total 2ml Vial Capacity: 630 Liquid Nitrogen Capacity: 34.8 Liters ½ CC Straw (1 Level Bulk): 3300 ½ CC Straw (10 Cane): 2100 Static Evaporation Rate (Liters/day): 0.18 Normal Working Time: 123 days Neck Diameter (mm): 88.9 Exterior Height (mm): 675.6 Exterior Diameter (mm): 462.3 Canister Diameter (mm): 71.1 Canister Height (mm): 279.4 Durable aluminum tank with powder coated exterior Lid lock provision fiberglass restricted opening neck design Compatible design for both liquid and vapor phase storage LN2 measuring rod for level validation Color-coded canisters for inventory identification One-year parts and labor warranty Three-year vacuum warranty </p>