

MFSD10 Polyclonal Antibody

Catalog number: 11518-1-AP


Size: 35 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

MFSD10; MFSD10, TETRA,
TETTRAN



mouse heart tissue were subjected to SDS PAGE followed by western blot with 11518-1-AP(MFSD10 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated MFSD10 MW:	455aa,48 kDa
Observed MFSD10 MW:	48-52 kDa
Positive WB detected in	Mouse heart tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2080
GenBank accession number:	BC014979
Gene ID (NCBI):	10227
Full name:	Major facilitator superfamily domain containing 10

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.