

# Nalgene Rapid-Flow Filters, Bottle Tops, and Analytical Filters

Get 20% off\* until year end (12/31/22)

## Nalgene Rapid-Flow units

Thermo Scientific™ Nalgene™ Rapid-Flow™ units feature the exclusive Rapid-Flow support column design, enabling fast flow rates and high throughput, thereby increasing filtration efficiency. They are an excellent choice for the filtration of cell culture media, sera, additives, and buffers.

### Polyethersulfone (PES) membrane

- The most common membrane material
- Fast flow rates and low clogging rates
- Available in 0.1 µm pore size for mycoplasma removal
- Tested for use with stem cells—a safe solution for maintaining sterility of stem cell cultures



## Nalgene Analytical Filters CN membranes

Use these economical sterile disposable filter funnels to perform microbiological QC testing and analysis of water, food/beverage, raw material and finished product. The squeeze and twist funnel is easy to remove for retrieval of the CN membrane.

- High impact polystyrene collars are rugged for reliable performance
- Pre-assembled, sterile, individually bagged to ensure sample integrity
- Adapters accommodate No. 8 stoppers for use with 0.5 in. hole in vacuum manifolds or filter flasks



## Ordering information

Product	Thomas No.	Mfr. No.	Promo Price
Nalgene Single Use Analytical Filter Funnels, CN, 0.2 µm, 100 mL, white, case of 50	4618U76	145-0020	\$233.60
Nalgene Single Use Analytical Filter Funnels, CN, 0.45 µm, 100 mL, white with black grid, case of 50	4618U7	145-0045	\$233.60
Nalgene Single Use Analytical Filter Funnels, CN, 0.2 µm, 250 mL, white, case of 50	4618U50	145-2020	\$260.00
Nalgene Single Use Analytical Filter Funnels, CN, 0.45 µm, 250 mL, white with black grid, case of 50	4618U53	145-2045	\$260.00
Nalgene Single Use Analytical Filter Funnels, CN, 0.45 µm, 100 mL, white with black grid, case of 50	4618U80	147-0045	\$245.60
Nalgene Single Use Analytical Filter Funnels, CN, 0.45 µm, 250 mL, white with black grid, case of 50	---	147-2045	---
Nalgene Sterile Analytical Filter Units, CN, 0.2 µm, 150 mL, pack of 12	1208M78	130-4020	\$555.20
Nalgene Sterile Analytical Filter Units, CN, 0.45 µm, 150 mL, case of 72	1208M79	130-4045	\$555.20
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.2 µm	1233Y81	564-0020	\$85.60
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.1 µm	1228Y92	565-0010	\$86.40
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.2 µm	1234K57	565-0020	\$88.00
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.45 µm	4618L93	165-0045	\$88.80
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.1 µm	1228Y93	568-0010	\$112.00
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.2 µm	1234K60	568-0020	\$113.60
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.45 µm	4618L97	168-0045	\$114.40
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.1 µm	1228Y94	566-0010	\$188.00
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.2 µm, 75 mm	1234K58	566-0020	\$192.00
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.2 µm, 90 mm	1234K61	569-0020	\$192.00
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.45 µm, 75 mm	4618M01	166-0045	\$192.00
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.45 µm, 90 mm	4618M03	169-0045	\$192.00
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.1 µm	1218X06	567-0010	\$254.40
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.2 µm	1234K59	567-0020	\$258.40
Nalgene Rapid-Flow Sterile Single Use Vacuum Filter Units, PES, 0.45 µm	4618M05	167-0045	\$260.80
Nalgene Rapid-Flow Sterile Single Use Bottle Top Filters A-PES, 45 mm, 500 mL, 75 mm, case of 12	1234K63	595-4520	\$116.80
Nalgene Rapid-Flow Sterile Single Use Bottle Top Filters A-PES, 45 mm, 150 mL, 50 mm	1234K65	596-4520	\$84.00
Nalgene Rapid-Flow Sterile Single Use Bottle Top Filters A-PES, 45 mm, 1 L, 90 mm	1234K67	597-4520	\$186.40
Nalgene Rapid-Flow Sterile Single Use Bottle Top Filters PES, 45 mm, 500 mL, 75 mm, 45 µm	4619N35	295-4545	\$116.80
Nalgene Rapid-Flow Sterile Single Use Bottle Top Filters PES, 45 mm, 150 mL, 50 mm, 45 µm	---	296-4545	---

Contact your Thomas Scientific sales representative for more information:

Name:

Email:

Phone:

**Get 20% off\* Rapid-Flow filters, bottle tops, and analytical filters**

\* Offer expires December 31, 2022. No additional notice of applied discounts will be given after checkout.

**thermo**scientific

**This product is intended for General Laboratory Use. It is the customer's responsibility to ensure that the performance of the product is suitable for customer's specific use or application.** © 2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **COL118799 0922**