

RFX5 Polyclonal Antibody

Catalog number: 12137-1-AP

Size: 28 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

RFX5; DNA binding protein

RFX5, Regulatory factor X 5,

RFX5



Jurkat cells were subjected to SDS PAGE followed by western blot with 12137-1-AP (RFX5 antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated RFX5 MW:	616aa, 65 kDa
Observed RFX5 MW:	65 kDa
Positive WB detected in:	Jurkat cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2783
GenBank accession number:	BC017471
Gene ID (NCBI):	5993
Full name:	Regulatory factor X, 5 (influences HLA class II expression)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.