



TheSharpChoice®

NMSDC Certified Diversity Supplier

DATASHEET

RELI® PRO Sutures

RELI® PRO Sutures from MYCO Medical are premium sterile absorbable and non-absorbable surgical-grade sutures packaged in a new low-memory racetrack tray to improve suture performance.

- All RELI® PRO sutures include needles made from 300 Series stainless steel for the perfect blend of strength and ductility
- RediPass™ needle coating facilitates smoother needle penetration and reduces tissue dragging
- 100% in-line manufacturing inspection ensures a smooth and secure needle to suture attachment site
- Packaged as 12 units per box, 30 boxes per case



Pure Glide

Continuous inner lined barrier eliminates any chance for the suture to get caught in gaps

Notched Label

Medical-grade paper label is notched in place to provide easy product identification and ensurable product integrity

Easy Lift Needle

The tab, when pressed from behind lifts the needle away from the card for easy removal

Twin Needle Parks

Securely accommodates and protects all MYCO needle configurations

RELI® PRO Suture Profiles


Absorbable Sutures



REDISORB® FAST PRO

Fast absorbing, coated, braided, Polyglycolic Acid (PGA), synthetic suture

- Undyed
- Similar to Ethicon VICRYL RAPIDE™



REDICAPRONE® 25 PRO

Monofilament, Poly(glycolide-co-caprolactone) (PGCL) synthetic suture

- Violet or Undyed
- Similar to Ethicon MONOCRYL®



REDISORB® PRO

Coated, braided, Polyglycolic Acid (PGA), synthetic suture


- Violet or Undyed
- Similar to Ethicon VICRYL®



REDIGUT® CHROMIC PRO

Chromic treated, twisted natural gut suture

- Natural
- Similar to Ethicon Chromic Gut



REDIDIOX® PRO

Monofilament, Polydioxanone (PDO), synthetic suture

- Violet
- Similar to Ethicon PDS®II



REDIGUT® PLAIN PRO

Plain, twisted, natural gut suture

- Natural
- Similar to Ethicon Plain Gut

Non-Absorbable Sutures



REDILON® PRO

Monofilament, nylon synthetic suture


- Black
- Similar to Ethicon ETHILON™



REDISILK® PRO

Organic silk natural suture

- Black
- Similar to Ethicon PERMAHAND® Silk



REDILENE® PRO

Monofilament, polypropylene, synthetic suture

- Blue
- Similar to Ethicon PROLENE™



REDIBOND® PRO

Coated, braided, polyester, synthetic suture

- Green
- Similar to Ethicon ETHIBOND EXCEL®

RELI® PRO Absorption Rates

Line	Material	Structure / Color	Strength Retention Profile	Time to Complete Absorption
REDISORB® FAST PRO	Fast-absorbing, synthetic Polyglycolic Acid (PGA)	Braided and Coated — Undyed	50% @ 5 days 0% @ 10–14 days	42 days
REDISORB® PRO	Synthetic Polyglycolic Acid (PGA)	Braided and Coated — Violet and Undyed	75% @ 14 days 50% @ 21 days	56–70 days
REDIDIOX® PRO	Synthetic Polydioxanone	Monofilament — Violet (dyed)	70% @ 14 days 60–50% @ 28 days	180–220 days
REDICAPRONE® 25 PRO	Synthetic Poly(glycolide-co-caprolactone) (PGCL)	Monofilament — Violet Monofilament — Undyed	80–60% @ 7 days 40–30% @ 14 days 60–50% @ 7 days 30–20% @ 14 days	90–120 days
REDIGUT® CHROMIC PRO	Chromic-treated Natural Gut	Twisted — Natural	Variable	Variable
REDIGUT® PLAIN PRO	Natural Gut	Twisted — Natural	Variable	Variable

The RELI® PRO Suture product catalog, with a full product code listings, and a complete brand cross reference chart are available upon request.

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com



MYCO
MEDICAL
TheSharpChoice®

2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com

www.mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2019 MYCO Medical, Inc. All rights reserved.

LIT-00593 Rev 0 Eff 09.23.2019