

POLR1D Polyclonal Antibody

Catalog number: 12254-1-AP

Size: 48 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

POLR1D; POLR1D, RPO13,

RPA9, POLR1C, hRPA19



HeLa cells were subjected to SDS PAGE followed by western blot with 12254-1-AP (POLR1D antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated POLR1D MW:	16 kDa
Observed POLR1D MW:	20kd
Positive WB detected in	HeLa cells, Jurkat cells, mouse lung tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2904
GenBank accession number:	BC015319
Gene ID (NCBI):	51082
Full name:	Polymerase (RNA) I polypeptide D, 16kDa

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.