

Uteroglobin Polyclonal Antibody

Catalog number: 10490-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

SCGB1A1; CC10, CC16, CCPBP,

CCSP, SCGB1A1, UGB, UP 1,

UP1, Urinary protein 1, Urine

protein 1, Uteroglobin

Applications

Tested applications:	ELISA
Cited applications:	FC, IF
Species specificity:	Human; other species not tested.
Cited species:	Human
Calculated Uteroglobin MW:	10 kDa
Observed Uteroglobin MW:	

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag0758
GenBank accession number:	BC004481
Gene ID (NCBI):	7356
Full name:	Secretoglobin, family 1A, member 1 (uteroglobin)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.