

## SF3B14 Polyclonal Antibody

Catalog number: 12379-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

SF3B14; CGI-110; HSPC175;

Ht006; P14; SAP14; SF3B14a



mouse uterus tissue were subjected to SDS PAGE followed by western blot with 12379-1-AP(SF3B14 Antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse uterus tissue
Calculated SF3B14 MW:	14kd
Observed SF3B14 MW:	14kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3045
GenBank accession number:	BC015463
Gene ID (NCBI):	51639
Full name:	Splicing factor 3B, 14 kDa subunit

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.