

## TIGD5 Polyclonal Antibody

Catalog number: 13644-1-AP


Size: 59 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

TIGD5; FLJ14926; MGC44883



HeLa cells were subjected to SDS PAGE followed by western blot with 13644-1-AP(TIGD5 Antibody) at dilution of 1:800

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HeLa cells
Calculated TIGD5 MW:	593aa,65kd
Observed TIGD5 MW:	70-80kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag4590
GenBank accession number:	BC032632
Gene ID (NCBI):	84948
Full name:	Tigger transposable element derived 5

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.