

## PDLIM3 Polyclonal Antibody

Catalog number: 13199-2-AP

Size: 61 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PDLIM3; ALP; DKFZp686L0362



human heart tissue were subjected to SDS PAGE followed by western blot with 13199-2-AP(PDLIM3 Antibody) at dilution of 1:1200

### Applications

|                          |  |
|--------------------------|--|
| Tested applications:     | ELISA, WB                                    |
| Species specificity:     | Human, mouse, rat; other species not tested. |
| Positive WB detected in: | Human heart tissue                           |
| Calculated PDLIM3 MW:    | 316aa,34kd                                   |
| Observed PDLIM3 MW:      | 39kd   |
| Recommended dilution:    | WB: 1:500-1:5000                             |

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

|                           |                      |
|---------------------------|----------------------|
| Immunogen:                | Ag3934               |
| GenBank accession number: | BC027870             |
| Gene ID (NCBI):           | 27295                |
| Full name:                | PDZ and LIM domain 3 |

### Product information

|                      |  |
|----------------------|--|
| Purification method: | Antigen affinity purification  |
| Storage:             | PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C. |