



Recommended Uses



General Description

Material	Vinyl	Type	Non-Sterile
Color	Clear	Usage	Single Use
Shape	Ambidextrous	Interior	Powder Free
Sizes	S - XL	Exterior	Smooth
Cuff	Beaded	Process	Polymer Coated

Product Codes

Small	20A00U981
Medium	20A00U982
Large	20A00U983
X-Large	20A00U984

Packaging

Box	100 Gloves Per Box
Box Dims	220 x 122 x 55 mm
Case	10 Boxes Per Case
Case Dims	293 x 254 x 220 mm

Physical Properties

Tensile Strength	11 Mpa (min)
Tensile Strength (Aged)	11 Mpa (min)
Powder Content	< 2 mg
FDA Listing	510(k) Status
Elongation	300%
Elongation (Aged)	300%
Protein Content	N/A
AQL Level	1.5

Dimensions

Weight	5.0 ± 0.2 g	Length	240 ± 5 mm
Thickness		Width	
Cuff	3 ± 1 mil	Small	85 ± 5 mm
Palm	3 ± 1 mil	Medium	95 ± 5 mm
Finger	4 ± 1 mil	Large	105 ± 5 mm
		X-Large	115 ± 5 mm

Barcodes

 6 97383 40101 4 S Box	 6 97383 40102 1 M Box	 6 97383 40103 8 L Box	 6 97383 40104 5 XL Box
 6 97383 41101 3 S Case	 6 97383 41102 0 M Case	 6 97383 41103 7 L Case	 6 97383 41104 4 XL Case