

ZEROTIP® Pipette Micro Tips

JET ZEROTIP® are manufactured from super clear high-quality Polypropylene. The surfaces of the tips are produced through a special process. This process makes the tip inner surface become super hydrophobic, thus significantly reduces sample loss and provides substantially higher reproducibility when working with critical media.



- Available with 6 transfer volumes of 10 µl, 20 µl, 100 µl, 200 µl, 300 µl and 1000 µl
- Free of Nucleic Acids, Pyrogens/Endotoxins, PCR Inhibitors and Trace Metals
- Ideal for biological samples that contain detergents such as SDS, Tween Triton® -100 ect
- Tips with PP filter are also available
- Preferred accessories for most brand micropipettor
- Non-Pyrogenic & DNase/RNase-free
- Packaging: 96/rack, 10 racks/box, 2 boxes/case



Thomas No.	Mfr. No.	Description	Pack/Case Qty.
21A00A321	PMT252010	ZEROTIP Universal pipette tip Sterile, Filter, 10 µl	960 pack/1920 Case
21A00A322	PMT233010	ZEROTIP Universal pipette tips, Sterile, Filter, 10XL µl	960 pack/1920 Case
21A00A323	PMT252020	ZEROTIP Universal pipette tips, Sterile, Filter, 20 µl	960 pack/1920 Case
21A00A324	PMT252100	ZEROTIP Universal pipette tips, Sterile, Filter, 100 µl	960 pack/1920 Case
21A00A325	PMT231200	ZEROTIP Universal pipette tips, Sterile, Filter, 200 µl	960 pack/1920 Case
21A00A326	PMT233300	ZEROTIP Universal pipette tips, Sterile, Filter, 300 µl	960 pack/1920 Case
21A00A327	PMT252000	ZEROTIP Universal pipette tips, Sterile, Filter, 1000 µl	960 pack/1920 Case
21A00A328	PMT371000	ZEROTIP Universal pipette tips, Sterile, Filter, 1000XL µl	960 pack/1920 Case