

## ARL4 Polyclonal Antibody

Catalog number: 10312-1-AP


Size: 30 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

ARL4A; ARL4, ARL4A



mouse liver tissue were subjected to SDS PAGE followed by western blot with 10312-1-AP (ARL4A antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated ARL4 MW:	23kd
Observed ARL4 MW:	24 kDa
Positive WB detected in	Mouse liver tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag0342
GenBank accession number:	BC001111
Gene ID (NCBI):	10124
Full name:	ADP-ribosylation factor-like 4A

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.