

POU2AF1 Polyclonal Antibody

Catalog number: 13605-1-AP


Size: 88 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

POU2AF1; POU2AF1, BOB1,
BOB 1, OBF1, OCAB



mouse lung tissue were subjected to SDS PAGE followed by western blot with 13605-1-AP(POU2AF1 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated POU2AF1 MW:	256aa,27 kDa
Observed POU2AF1 MW:	35kd
Positive WB detected in	Mouse lung tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4482
GenBank accession number:	BC032549
Gene ID (NCBI):	5450
Full name:	POU class 2 associating factor 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.