

SAMSN1 Polyclonal Antibody

Catalog number: 13063-1-AP


Size: 124 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

SAMSN1; HACS1, NASH1,
SAMSN1, SASH2, SH3 SAM
adaptor protein, SH3D6B



Immunofluorescent analysis of HepG2 cells using 13063-1-AP (SAMSN1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



K-562 cells were subjected to SDS PAGE followed by western blot with 13063-1-AP (SAMSN1 Antibody) at dilution of 1:1500

Applications

Tested applications:	ELISA, WB, IF
Cited applications:	IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Cited species:	Human
Calculated SAMSN1 MW:	42 kDa
Observed SAMSN1 MW:	42-50 kDa
Positive WB detected in	K-562 cells
Positive IF detected in	HepG2 cells
Recommended dilution:	WB: 1:500-1:5000 IF: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3755
GenBank accession number:	BC029112
Gene ID (NCBI):	64092
Full name:	SAM domain, SH3 domain and nuclear localization signals 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.