

GSTA2 Polyclonal Antibody

Catalog number: 10505-1-AP


Size: 60 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

GSTA2; Glutathione S transferase A2, GST class alpha member 2, GST gamma, GST HA subunit 2, GST2, GSTA2, GSTA2 2, GTA2, GTH2



Immunohistochemical of paraffin-embedded human pancreas cancer using 10505-1-AP (GSTA2 antibody) at dilution of 1:100 (under 10x lens)



HEK-293 cells were subjected to SDS PAGE followed by western blot with 10505-1-AP (GSTA2 antibody) at dilution of 1:100

Applications

Tested applications:	ELISA, WB, IHC
Cited applications:	WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Cited species:	Mouse
Calculated GSTA2 MW:	26 kDa
Observed GSTA2 MW:	26kd
Positive WB detected in:	HEK-293 cells, HeLa cells
Positive IHC detected in:	Human pancreas cancer tissue
Recommended dilution:	WB: 1:50-1:200 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag0781
GenBank accession number:	BC002895
Gene ID (NCBI):	2939
Full name:	Glutathione S-transferase alpha 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.