

## AMPK beta 2 Polyclonal Antibody

Catalog number: 14429-1-AP


Size: 29 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PRKAB2; AMPK beta 2, AMPK subunit beta 2, PRKAB2



human brain tissue were subjected to SDS PAGE followed by western blot with 14429-1-AP (PRKAB2 antibody) at dilution of 1:400

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated AMPK beta 2 MW:	30 kDa
Observed AMPK beta 2 MW:	34 kDa
Positive WB detected in	Human brain tissue, human heart tissue, mouse testis tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag5806
GenBank accession number:	BC053610
Gene ID (NCBI):	5565
Full name:	Protein kinase, AMP-activated, beta 2 non-catalytic subunit

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.