

GUK1 Polyclonal Antibody

Catalog number: 10755-1-AP


Size: 38 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

GUK1; GMK, GMP kinase,
Guanylate kinase, guanylate
kinase 1, GUK1



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with 10755-1-AP (GUK1 Antibody) at dilution of 1:350

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated GUK1 MW:	22 kDa
Observed GUK1 MW:	
Positive WB detected in	Transfected HEK-293 cells
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1205
GenBank accession number:	BC009914
Gene ID (NCBI):	2987
Full name:	Guanylate kinase 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.