

## UL27 Polyclonal Antibody

Catalog number: 10363-1-AP

Size: 150 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

;

### Applications

**Tested applications:** ELISA  
**Species specificity:** Hsv; other species not tested.  
**Calculated UL27 MW:** 69kd  
**Observed UL27 MW:**

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

**Immunogen:** Ag0411  
**GenBank accession number:** NP\_044629  
**Gene ID (NCBI):**  
**Full name:** Undefined

### Product information

**Purification method:** Antigen affinity purification  
**Storage:** PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C, avoid freeze/thaw cycles.