

MXI1 Polyclonal Antibody

Catalog number: 12360-1-AP

Size: 21 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

MXI1; bHLHc11, MAD2, Max

interacting protein 1, MAX

interactor 1, MXD2, MXI, MXI1



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 12360-1-AP (MXI1 Antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated MXI1 MW:	419aa, 47 kDa
Observed MXI1 MW:	29 kDa
Positive WB detected in:	SH-SY5Y cells, mouse brain tissue
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3024
GenBank accession number:	BC012907
Gene ID (NCBI):	4601
Full name:	MAX interactor 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.