


GEB Universal tips for TWC thin wallcrown And for improved flexibility.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (1.8*3.5)mm	Packaging Spec.
PT0010-B-N	10ul	Natural				1000tips/bag, 20bags/case, 20000tips/case
PT0010-B-N-LB	10ul	Natural	Low Retention			
Case dimension L*W*H: 455*265*255mm. Gross Weight: 7.6kg.						
FT0010-B-N	10ul	Natural			Filtered	1000tips/bag,
FT0010-B-N-LB	10ul	Natural	Low Retention		Filtered	10bags/case,
FT0010-B-NS	10ul	Natural		Sterilized	Filtered	10000tips/case
FT0010-B-NS-LB	10ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 450*270*190 mm. Gross Weight: 4.5kg.						
PT0010-9-N	10ul	Natural				96tips/rack,
PT0010-9-N-LB	10ul	Natural	Low Retention			10racks/pack,
PT0010-9-NS	10ul	Natural		Sterilized		5packs/case,
PT0010-9-NS-LB	10ul	Natural	Low Retention	Sterilized		4800tips/case.
FT0010-R-NS	10ul	Natural		Sterilized	Filtered	
FT0010-R-NS-LB	10ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 480*345*280mm. Gross Weight: 7.6kg.						

- ✧ 10μL with graduation marks at 2ul and 10ul
- ✧ 96 Tip/Hanging Tip Rack format / Bulk
- ✧ Polypropylene
- ✧ Presterile/Nonsterile
- ✧ Filtered (With aerosol-resistant filters) /NonFiltered
- ✧ Certified RNase-/DNase-free
- ✧ Nonpyrogenic
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Universal tips workstations.

Drawing:



Prepared by: Anne Oct.19, 2020 Approve: WeiMing Yu Oct.19,2020