

IL37 Polyclonal Antibody

Catalog number: 11863-1-AP


Size: 27 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

IL1F7; FIL1, FIL1 zeta,
FIL1(ZETA), FIL1Z, IL 1 zeta, IL
1F7, IL 1H4, IL 1RP1, IL 1X,
IL1F7, IL1H4, IL1RP1,
Interleukin 1 family member 7,
Interleukin 1 homolog 4,
Interleukin 1 related protein,
Interleukin 1 zeta



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with 11863-1-AP(IL1F7 Antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated IL37 MW:	218aa,24 kDa
Observed IL37 MW:	
Positive WB detected in	Transfected HEK-293 cells, mouse lung tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2441
GenBank accession number:	BC020637
Gene ID (NCBI):	27178
Full name:	Interleukin 1 family, member 7 (zeta)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.