

SPANXC Polyclonal Antibody

Catalog number: 11957-1-AP


Size: 75 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

SPANXC; C, Cancer/testis antigen 11.3, CT11.3, CTp11, SPANX C, SPANX family member C, SPANX family, member C, SPANXC



human testis tissue were subjected to SDS PAGE followed by western blot with 11957-1-AP (SPANXC antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated SPANXC MW:	97aa, 11 kDa
Observed SPANXC MW:	20 kDa
Positive WB detected in:	Human testis tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2561
GenBank accession number:	BC054023
Gene ID (NCBI):	64663
Full name:	SPANX family, member C

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.