

TECHNICAL DATA SHEET


Meets
USP <800>
Guidelines

Tested
against
Fentanyl,
Xylazine
& Gastric
Acid

Silicone-
Free &
Motor Oil
Resistant

Description

HALYARD* **PUREZERO*** BLACK (9.5") Nitrile Exam Gloves are designed for Medical, Life Sciences, Forensics, Automotive, Janitorial, Food & Beverage, Industrial, Machinery & Equipment, EMS, Utilities, Warehousing, Healthcare and Laboratories. HALYARD* **PUREZERO*** BLACK Nitrile Exam Gloves are tested for use against 15 chemotherapy drugs and 45 chemicals. These gloves have also been tested and shown to have low dermatitis potential and no added/non-detectable levels of residual chemicals, including accelerators and contact sensitizers. The beaded cuff allows for ease of donning, and fully textured glove surface provide tactile grip.




Physical Properties

AQL	1.5
Non-Sterile	✓
Ambidextrous	✓
Fully Textured	✓
Not Made With Natural Rubber Latex	✓
Powder-Free	✓
Tensile Strength	18 MPa (Target)
Ultimate Elongation	500%
Shelf Life	3 Years

SMARTPULL* Dispenser

PUREZERO* BLACK gloves have a SMARTPULL* Dispenser, which incorporates two separate openings on the box. The first, smaller opening is used when the box is full to control dispensing and lessen waste. When the box is half empty, the second, larger opening allows easier access to the gloves.



Glove Dimensions

	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE
BLACK Glove Length (inch/cm)	9.5"/24	9.5"/24	9.5"/24	9.5"/24	9.5"/24
Middle Finger Length (mm)	77	78	85	88	88
Width of Palm (mm)	70	80	95	106	116
Fingertip Thickness	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)	0.13 mm (5.1 mil)
Palm Thickness	0.11 mm (4.3 mil)	0.11 mm (4.3 mil)	0.11 mm (4.3 mil)	0.11 mm (4.3 mil)	0.11 mm (4.3 mil)
Cuff Thickness	0.09 mm (3.5 mil)	0.09 mm (3.5 mil)	0.09 mm (3.5 mil)	0.09 mm (3.5 mil)	0.09 mm (3.5 mil)



HALYARD* PUREZERO* BLACK Nitrile Exam Gloves For Medical/Industrial/Lab Use

TECHNICAL DATA SHEET

Ordering Information

Product Name	Thomas #	Mfgr #	Size	Box Qty.	Case Qty.
HALYARD* PUREZERO* BLACK Nitrile Exam Gloves - 9.5"/24cm Length	25A00A970	LFS711XS	XS	100	1000
	25A00A971	LFS711SM	S	100	1000
	25A00A972	LFS711MD	M	100	1000
	25A00A973	LFS711LG	L	100	1000
	25A00A974	LFS711XL	XL	90	900

Storage Instructions

HALYARD* PUREZERO* Nitrile Gloves shall be stored in conditions where the product is kept dry (away from moisture), away from direct sunlight, away from sources of heat, and away from radioactive sources.

Quality & Regulatory Standards

Compliant to quality management systems:

- ISO 9001: 2015
- ISO 13485: 2016
- US FDA 21 CFR part 820

Compliant to these biocompatibility standards:

- ISO 10993-10
- ISO 10993-11

Compliant to these food handling regulations:

- US FDA 21 CFR 177.2600 and US FDA CFR 180.22

Compliant to these medical device regulations:

- US FDA 21 CFR 880.6250 (510 (k) cleared device)

Compliant to these standards:

- ASTM D6319-19:2023
- ASTM F739-20:2020
- ASTM D6978-05: 2019
- ASTM F1671-07:2010
- Modified Draize - 95
- ASTM E1252-98: 2021



LOW
DERMATITIS
POTENTIAL



For additional information
or samples, contact your
local distributor.

This fact sheet has been created using the most recent information. In the interest of continuous improvement, the characteristics of the products may change without prior notice.

*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2025. All rights reserved.
009099 COPY-05769 November 2025



HALYARD* is an Owens & Minor brand



customerservice@thomassci.com



thomassci.com

