

# KingFisher™ Compatible Tip Combs

- Low-binding, medical grade polypropylene
- Compatible with KingFisher™ & Accuris IsoPure™
- 8, 24, and 96 tip configurations
- Available sterile or non-sterile

labForce® tip combs are produced from low-binding, medical grade polypropylene and are made specifically to fit KingFisher™ Flex, Accuris IsoPure™ and a variety of other automation platforms.

Applications include DNA extraction, quantification, viral RNA extraction, immunoprecipitations, etc.



## KingFisher™ Compatible Tip Combs - *Sterile*

Cat. No.	Description	Qty
23A00W741	Tip Comb, 8-Tip, V-Bottom, PP, for KingFisher™, Sterile, 240 Bags of 2	480
23A00W743	Tip Comb, 24-Tip, V-Bottom, PP, for KingFisher™, Sterile, 25 Bags of 2	50
23A00W745	Tip Comb, 96-Tip, V-Bottom, PP, for KingFisher™, Sterile, 50 Bags of 2	100

## KingFisher™ Compatible Tip Combs - *Non-Sterile*

Cat. No.	Description	Qty
23A00W740	Tip Comb, 8-Tip, V-Bottom, PP, for KingFisher™, 240 Bags of 2	480
23A00W742	Tip Comb, 24-Tip, V-Bottom, PP, for KingFisher™, 25 Bags of 2	50
23A00W744	Tip Comb, 96-Tip, V-Bottom, PP, for KingFisher™, 50 Bags of 2	100

\*KingFisher™ is a trademark of Thermo Fisher Scientific