

**BIOFIL®**

Document Number

Revision

Page

Date

**GSP002****A0**

1 of 1

01-24-2018

**PRODUCT DATA SHEET****Serological Pipets****2mL Serological Pipets**

Description: 2mL in 1/50mL


Purpose: Use for liquid transfer

**Materials**

Pipet: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Color: Clear



Unit: mm

	Dimensions	Tolerance
A	∅6.3	±0.1
B	∅2.93	±0.35
C	272.0	±2.0

**Features**

- Temperature: -20°~+50°
- Sterilization: sterilized by gamma irradiation SAL10<sup>6</sup>(ISO11137)
- Latex Statement: this product is latex free
- DNase/RNase free, non-pyrogenic, non-cytotoxic
- Meet the requirement of BSE/TSE
- Use dispenser box for easy opening
- Available in Individually packed in peel-to-open wrap
- Autoclavability: No
- Negative graduations up to -0.6 ml (total volume = 2.6 ml)
- Graduation (mL): in 1/50mL
- Graduations are calibrated for accurate dispensing to within ±2%
- Color coded: green
- Every product is printed with a batch number for traceability
- Shelf Life: 3 years after month of production
- Manufactured in a class 100,000 room environment
- Manufactured under EN ISO13485:2003 and ISO9001:2008 quality management system

②

**2mL Serological Pipets**

Thomas No.	Cat.No	Description	Sterile	Qty/Tray	Qty/Pk	Qty/Case
21A00B814	GSP010002	2mL in 1/50mL, Individually Wrapped, Paper/Plastic	Yes	1	500	500

