

## ICP27 Polyclonal Antibody

Catalog number: 10186-1-AP

Size: 60 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

ICP27, Immediate-early protein

IE63, Infected cell protein 27,

mRNA export factor, UL54,

VMW63

### Applications

**Tested applications:** ELISA  
**Species specificity:** HSV; other species not tested.  
**Calculated ICP27 MW:** 55 kDa  
**Observed ICP27 MW:**

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

**Immunogen:** Ag0239  
**GenBank accession number:** AAF43147  
**Gene ID (NCBI):**  
**Full name:** Undefined

### Product information

**Purification method:** Antigen affinity purification  
**Storage:** PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.