

Product Name:	ABS 1
Description:	<p> ABS CryoMax System Vial Capacity: 750 (2ml) Total Box Capacity: 30 Rack Capacity: 6 Boxes per rack: 5 Vials per box: 25 Liquid Nitrogen Capacity: 47.4 liters Static Evaporation Rate: .39 liters/day Static Holding Time: 76 days Neck Diameter (mm): 127 Exterior Height (mm): 673.1 Exterior Diameter (mm): 508 Durable aluminum tank with powder coated exterior Stainless steel storage racks, configured for 2" high boxes Roller base Liquid nitrogen measuring rod for level validation Lid lock provision Insulated neck tube core Compatible design for both liquid and vapor phase storage One-year warranty on all parts and labor Three-year vacuum warranty </p>