

## TFAP2B Polyclonal Antibody

Catalog number: 13183-1-AP


Size: 27 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

TFAP2B; AP 2B, AP2 B, AP2 beta, TFAP2B, Transcription factor AP 2 beta



MCF7 cells were subjected to SDS PAGE followed by western blot with 13183-1-AP(TFAP2B antibody) at dilution of 1:800



IP Result of anti-TFAP2B (IP:13183-1-AP, 4µg; Detection:13183-1-AP 1:500) with MCF-7 cells lysate 800µg.

### Applications

Tested applications:	ELISA, WB, IP
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated TFAP2B MW:	449aa, 49 kDa
Observed TFAP2B MW:	50kd
Positive WB detected in	MCF7 cells
Positive IP detected in	MCF-7 cells
Recommended dilution:	WB: 1:200-1:2000 IP: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3853
GenBank accession number:	BC037225
Gene ID (NCBI):	7021
Full name:	Transcription factor AP-2 beta (activating enhancer binding protein 2 beta)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.