

STK24 Polyclonal Antibody

Catalog number: 13579-1-AP

Size: 61 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

STK24; MST-3; MST3; MST3B;

STE20; STK3



HeLa cells were subjected to SDS PAGE followed by western blot with 13579-1-AP(STK24 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HeLa cells
Calculated STK24 MW:	431aa,48kd
Observed STK24 MW:	52kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4483
GenBank accession number:	BC035578
Gene ID (NCBI):	8428
Full name:	Serine/threonine kinase 24 (STE20 homolog, yeast)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.