



# ACCELERATOR SULFUR FREE NITRILE CLASS 100 ISO 5 12" GLOVES

## NCP LITE SERIES

Maximize cleanroom safety with our CTI NCP Lite Series Class 100 12" Nitrile Gloves. Accelerator- and sulfur-free, these gloves reduce contamination risks and allergic reactions, providing durable, comfortable protection for critical environments.



<b>100% Nitrile - Accelerator and Sulfur Free</b>	
Class 100 (ISO 5) Clean Room Compatible or Higher	
No Latex Proteins	
White Color, Super Soft Feel	
Anti-Static	
Textured for Critical Handling	
Double Bagged	
Tapered, Beaded Cuff	
12" Length, Ambidextrous	
Leak Tested	
<b>Length:</b>	12" (+/-0.5)
<b>Thickness:</b>	3.5 mil (+/-2mm)
<b>Tensile Strength:</b>	2175 psi/15 mpa
<b>Elongation:</b>	600%
<b>Packaging:</b>	Poly Bagged 100/Pack, 1000/Case
<b>Thomas No. (Cat No.)</b>	24A00G489 (NCPL201)
<b>Sizes:</b>	Small
	24A00G490 (NCPL202)
	Medium
	24A00G491 (NCPL203)
	Large
	24A00G492 (NCPL204)
	X-Large

Meets or Exceeds ASTM D3578-95 Standards  
Particle and Extractable Testing in Accordance with IES-RP-CC-005.2  
Sampled in Accordance with MIL STD-105E  
Qualified for Use in the Pharmaceutical Industry



ThomasSci.com  
833.544.SHIP (7447)  
CustomerService@thomassci.com

Connect With Us:

