

## OXR1 Polyclonal Antibody

Catalog number: 13514-1-AP


Size: 44 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

OXR1; FLJ10125, FLJ38829,  
FLJ41673, FLJ42450, FLJ45656,  
Nbla00307, oxidation resistance  
1, Oxidation resistance protein  
1, OXR1



K-562 cells were subjected to SDS  
PAGE followed by western blot with  
13514-1-AP( OXR1 Antibody) at  
dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated OXR1 MW:	758aa, 85 kDa
Observed OXR1 MW:	120-140 kDa
Positive WB detected in:	K-562 cells, HeLa cells, Neuro-2a cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag4439
GenBank accession number:	BC032710
Gene ID (NCBI):	55074
Full name:	Oxidation resistance 1

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.