

- VersaCap™ technology provides the flexibility of interchangeable top-connectors
- VersaCap™ are offered in 50mm, 80mm, and 120mm industry standard sizes
- Quick connect adapters available for easy port changes
- Leakproof cap with replaceable gasket
- VersaCap™ are stackable
- Available in a variety of configurations
- VersaCap™, VersaBarb™, and Spigots are autoclavable
- VersaBarb™ offer a convenient solution for bottom-connections and can fit a variety of hose sizes
- Spigot offers two flow rates

Foxx Life Sciences has enhanced the features and benefits of the EZgrip™ carboys and Heavy Duty Vacuum Bottles by offering the VersaCap™, VersaBarb™ and state-of-the-art spigot technology. These three components expand the wide array of applications supported by the EZgrip™ carboy.

Patented VersaCap™ are leakproof, have replaceable platinum-cured silicone gaskets, and provide the ultimate in top-connection flexibility. Open VersaCap™ are compatible with a variety of adapter inserts, enabling users to quickly change the top connections of each carboy or bottle to match their application needs. VersaCap™ offer large handles with finger grips that allow for easy tightening and opening. They are also stackable, and are offered in 50 mm, 80 mm (compatible with most carboys) and 120 mm wide-mouth versions that fit exclusively on Roush Life Sciences' 20L and 40L EZfil™ carboys.

Quick-connect hose barb options include quick-connect port combinations of two, three, and four, with 1/4" and 3/8" hose barb sizes, as well as molded-in hose barb options, available with 1/8", 1/4", 5/16", and 1/2" sizes. The ClearView™ adapter is also available for users who need a crystal-clear window into the carboy for visualizing solution homogenization.

For bottom tubulation options, VersaBarb™ provide flexibility in drainage, pump connection, filling bell, and automated equipment. VersaBarb™ are available in 1/4", 3/8", 1/2", and 3/4" hose barb sizes. A 3/4" sanitary fitting VersaBarb™ is also available to enable the EZgrip™ carboy to function as a supply reservoir to a larger system.


Foxx Life Sciences' EZgrip™ carboy spigots have large levers that control two flow rates, including a continuous flow option. The spigot spout features 5/8" and 3/4" hose barbs, and a standard internal luer taper to fit a male barb adapter. The spigot body has internal threads with a hard stop to prevent over-tightening which may damage the seal.

All VersaCap™, adapters, VersaBarb™, and spigots are made from polypropylene, which may be autoclaved at 15 PSI and 121°C (250°F) for 20 minutes.








VersaCap™ 50mm


	205-0313-RLS			
	Cap, Closed			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	n/a	n/a	2

	205-0345-RLS			
	Cap, Vented top with 0.2 μm Hydrophobic Membrane			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	n/a	n/a	1


	205-0314-RLS			
	Cap, Open with closed adapter			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	n/a	n/a	2


	205-3331-RLS			
	Adapter, Closed			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	n/a	n/a	1


	205-0317-RLS			
	Cap, Open top with adapter, 1/8" HB x 2 and Vent			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	1/8"	n/a	1

	205-5331-RLS			
	Adapter, 1/8" HB x 2 and Vent			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	1/8"	n/a	1

	205-0318-RLS			
	Cap, Open top with adapter, 1/4" HB x 2 and Vent			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	1/4"	n/a	1

	205-5332-RLS			
	Adapter, 1/4" HB x 2 and Vent			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	1/4"	n/a	1

	205-0319-RLS			
	Cap, Open top with adapter, 5/16" HB x 2 and Vent			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	5/16"	n/a	1


	205-5333-RLS			
	Adapter, 5/16" HB x 2 and Vent			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	50	5/16"	n/a	1


HB = Hose Barb








VersaCap™ 80mm


	205-1003-RLS Cap, Closed	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	n/a	n/a	2


	205-1004-RLS Cap, Open with closed adapter	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	n/a	n/a	2


	205-1007-RLS Cap, Open top with adapter, 1/4" HB x 2 and Vent	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	1/4"	n/a	1


	205-1008-RLS Cap, Open top with adapter, 1/2" HB x 2 and Vent	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	1/2"	n/a	1


	205-3001-RLS Adapter, Closed	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	n/a	n/a	1


	205-3002-RLS Adapter, 1/4" HB x 2	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	1/4"	2	1


	205-3003-RLS Adapter, 3/8" HB x 2	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	3/8"	2	1


	205-3004-RLS Adapter, 1/4" HB x 3	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	1/4"	3	1


	205-3005-RLS Adapter, 3/8" HB x 3	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	3/8"	3	1

	205-3006-RLS Adapter, 1/4" HB x 4	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	1/4"	4	1

	205-3007-RLS Adapter, 3/8" HB x 4	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	3/8"	4	1

	205-5001-RLS Adapter, 1/4" HB x 2 and Vent	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	1/4"	n/a	1


	205-5002-RLS Adapter, 1/2" HB x 2 and Vent	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	1/2"	n/a	1


	208-3001-RLS Adapter, Closed with Clear Window	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
		80	n/a	n/a	1


HB = Hose Barb





VersaCap™ 120mm


	205-2003-RLS			
	Cap, Closed			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	n/a	n/a	2


	205-4004-RLS			
	Adapter, 1/4" HB x 3			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	1/4"	3	1


	205-2004-RLS			
	Cap, Open with closed adapter			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	n/a	n/a	2


	205-4005-RLS			
	Adapter, 3/8" HB x 3			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	3/8"	3	1


	205-4001-RLS			
	Adapter, Closed			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	n/a	n/a	1

	205-4006-RLS			
	Adapter, 1/4" HB x 4			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	1/4"	4	1

	205-4002-RLS			
	Adapter, 1/4" HB x 2			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	1/4"	2	1

	205-4007-RLS			
	Adapter, 3/8" HB x 4			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	3/8"	4	1

	205-4003-RLS			
	Adapter, 3/8" HB x 2			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	3/8"	2	1

	205-4008-RLS			
	Adapter, EZ-Top adapter			
	Cap Size	Hose Barb Size	No. of Quick Connects	Qty/pk
	120	n/a	n/a	1

HB = Hose Barb





VersaBarb™ and Spigot Technology


225-4001-RLS

VersaBarb™, Spigot fitting, closed cap

Hose Barb Size	Thread Diameter	Qty/pk
n/a	1 1/8"	1


225-4005-RLS

VersaBarb™, Spigot fitting, with 3/4" HB

Hose Barb Size	Thread Diameter	Qty/pk
3/4"	1 1/8"	1


225-4002-RLS

VersaBarb™, Spigot fitting, with 1/4" HB

Hose Barb Size	Thread Diameter	Qty/pk
1/4"	1 1/8"	1


225-4006-RLS

VersaBarb™, Spigot fitting, 3/4" Sanitary Connector

Hose Barb Size	Thread Diameter	Qty/pk
n/a	1 1/8"	1


225-4003-RLS

VersaBarb™, Spigot fitting, with 3/8" HB

Hose Barb Size	Thread Diameter	Qty/pk
3/8"	1 1/8"	1


225-1001-RLS

Spigot, 2 pack

Hose Barb Size	Thread Diameter	Qty/pk
n/a	1 1/8"	2


225-4004-RLS

VersaBarb™, Spigot fitting, with 1/2" HB

Hose Barb Size	Thread Diameter	Qty/pk
1/2"	1 1/8"	1

HB = Hose Barb



We recommend testing for your specific application before use.