



Specifications

Electrical requirements	208-240VAC 50/60Hz single phase, 3800 watts
Chamber dimensions	27.5 x 25 x 35.5" (WxDxH)
Unit dimensions	33.5 x 34 x 52" (WxDxH)
Shipping dimensions	40 x 38 x 58" (WxDxH)
Temperature range	Ambient temperature to 300°C
Temperature uniformity	+/- 2°F
Temperature controller precision	0.1°F
Timer range	1 to 9999 minutes
Racks included	2 (additional racks available at extra cost)
Maximum rack capacity	16 (height between each rack is 2 inches)
Exhaust port	2" diameter x 2, adjustable
Shipping weight	435 Lb
Access port	2" diameter with cover, one
Certification	CE, UL/CSA (optional)

Standard package

Part description	Quantity	Part image
Ai 14.5 cu ft forced air oven	1 pc	
Stainless steel rack and brackets	2 sets	
User's manual	1 pc	