



Image shown for reference purposes only. Actual product appearance may vary.

SOFT TOUCH* II Surgical Mask, Horizontal Ties, Blue

47500

Halyard Health Face Masks offer the broadest range of styles, features and levels of protection available today. The materials and donning attachments are sonically bonded, and all masks have enclosed nosepieces to assist in conforming to the contours of the face. All face masks are natural rubber latex-free.

- Pleat style with horizontal ties
- Protective three-layer construction
- Appropriate for situations where exposure to blood and/or bodily fluids is not a risk.

Items Per Case	300
Boxes Per Case	6
Items Per Box	50

See Packaging Specifications

Packaging Details:	1 CS = 6 DSP	Cases per Layer:	8 Case(s) per Layer
	6 DSP = 50 Units	Cases per Stack:	14 Case(s) per Stack
	1 CS = 300 Units	Layers per Stack:	112 Layer(s) per Stack
Volume:	0.964 FT ³	Pallet Load Specifications:	GMA: 48.000" x 40.000" x 46.140"
Length:	18.440 IN		SGM: 48.000" x 40.000" x 23.070"
Width:	11.750 IN		TL: 48.000" x 40.000" x 99.970"
Height:	7.690 IN	Floor Load Specifications:	GMA: 47.000" x 36.880" x 53.830"
Weight:	5.080 LB		SGM: 47.000" x 36.880" x 30.760"

Specifications

Brand	HALYARD*
UNSPSC Classification	42131606
Sterile	No
Color	Blue
Name	SOFT TOUCH* II Surgical Mask, Horizontal Ties, Blue
UK SKU	P4750001
Type	Surgical
Style	Pleat
Fluid Resistant	No
Fluid Protection	None
BFE (3)	≥ 98%
SPLASH RESISTANCE (80,120,160)	0
PFE (.1)	≥ 97%
Items Per Box	50
Boxes Per Case	6
Items Per Case	300

Copyright 2018 Halyard Worldwide, Inc. All Rights Reserved. Your visit to this site and use of the information hereon is subject to the terms of our Legal Statement. Please Review our Privacy Statement. Please Review our California Compliance Declaration. *Registered Trademark or Trademark of Halyard Worldwide, Inc., or its affiliates.