

WHITE DIAMOND SERIES BUILT-IN UNDERCOUNTER REFRIGERATORS AND FREEZERS

Refrigerators

Model	Cu. Ft.	Defrost	Door	Shelves	W"	D"	H"	RL	H.P.	Amps
PR051WWW/0	4.6	Cycle	1 Solid	2 Adj./ 1 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.74
PR051WWWLH/0	4.6	Cycle	1 Solid LH	2 Adj./ 1 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.74
PR051WWG/0	4.6	Cycle	1 Glass	2 Adj./ 1 Fixed	23 3/4	25 1/2	33 3/8	R600a	1/8	1.74
PR051WWGLH/0	4.6	Cycle	1 Glass LH	2 Adj./ 1 Fixed	23 3/4	25 1/2	33 3/8	R600a	1/8	1.74
PR051WWGADA/0	4.6	Cycle	1 Glass	2 Adj./ 1 Fixed	23 3/4	25 1/2	31 15/16	R600a	1/8	1.74
PR051WWGADALH/0	4.6	Cycle	1 Glass LH	2 Adj./ 1 Fixed	23 3/4	25 1/2	31 15/16	R600a	1/8	1.74
PR051WWWADA/0	4.6	Cycle	1 Solid	2 Adj./ 1 Fixed	23 3/4	25 1/2	31 15/16	R600a	1/8	1.74
PR051WWWADALH/0	4.6	Cycle	1 Solid LH	2 Adj./ 1 Fixed	23 3/4	25 1/2	31 15/16	R600a	1/8	1.74
PR051SSS/0	4.5	Cycle	1 Solid	2 Adj./ 1 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.74
PR051SSSLH/0	4.5	Cycle	1 Solid LH	2 Adj./ 1 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.74

Freezers

Model	Cu. Ft.	Defrost	Door	Shelves	W"	D"	H"	RL	H.P.	Amps
PF051WWW/0M	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.5
PF051WWWLH/0M	4.2	Manual	1 Solid LH	2 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.5
PF051WWWADA/0M	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.5
PF051WWWADALH/0M	4.2	Manual	1 Solid LH	2 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.5
PF051SSS/0M	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	31 15/16	R600a	1/8	1.5
PF051SSSLH/0M	4.2	Manual	1 Solid LH	2 Fixed	23 3/4	24 1/2	31 15/16	R600a	1/8	1.5
PF051WWW/0A	4.2	Auto	1 Solid	3 Adj.	23 3/4	24 1/2	33 3/8	R600a	1/8	1.8
PF051WWWLH/0A	4.2	Auto	1 Solid LH	3 Adj.	23 3/4	24 1/2	33 3/8	R600a	1/8	1.8

WHITE DIAMOND SERIES FREESTANDING UNDERCOUNTER REFRIGERATORS AND FREEZERS

Refrigerators

Model	Cu. Ft.	Defrost	Door	Shelves	W"	D"	H"	RL	H.P.	Amps
PR011WWW/0	1	Cycle	1 Solid	2 Adj.	17 1/4	19 1/4	21 1/4	R600a	1/8	.90
PR011WWWLH/0	1	Cycle	1 Solid LH	2 Adj.	17 1/4	19 1/4	21 1/4	R600a	1/8	.90
PR011WWG/0	1	Cycle	1 Glass	2 Adj.	17 1/4	19 1/4	21 1/4	R600a	1/8	.90
PR011WWGLH/0	1	Cycle	1 Glass LH	2 Adj.	17 1/4	19 1/4	21 1/4	R600a	1/8	.90
PR031WWW/0	2.5	Cycle	1 Solid	4 Adj.	17 3/4	19 1/2	28 3/8	R600a	1/8	.90
PR031WWG/0	2.5	Cycle	1 Glass	4 Adj.	17 3/4	19 1/2	28 3/8	R600a	1/8	.90
PR061WWW/0	5.2	Cycle	1 Solid	3 Adj.	23 3/4	24	32 1/8	R600a	1/8	1.74
PR061WWG/0	5.2	Cycle	1 Glass	3 Adj.	23 3/4	24	32 1/8	R600a	1/8	1.74

Freezers

Model	Cu. Ft.	Defrost	Door	Shelves	W"	D"	H"	RL	H.P.	Amps
PF021WWW/0M	1.7	Manual	1 Solid	2 Adj.	23 3/4	24	20 3/8	R600a	1/8	1
PF021WWWLH/0M	1.7	Manual	1 Solid LH	2 Adj.	23 3/4	24	20 3/8	R600a	1/8	1
PF021WWW/0A	1.3	Auto	1 Solid	2 Adj.	23 3/4	24	20 3/8	R600a	1/8	1
PF021WWWLH/0A	1.3	Auto	1 Solid LH	2 Adj.	23 3/4	24	20 3/8	R600a	1/8	1
PF051WWW/0M30	4	Manual	1 Solid	3 Adj.	23 3/4	24	32 1/8	R600a	1/4	2.2
PF051WWW/0M40	4	Manual	1 Solid	3 Adj.	23 3/4	24	32 1/8	R600a	1/4	2.2

WHITE DIAMOND SERIES FLAMMABLE STORAGE UNDERCOUNTER REFRIGERATOR AND FREEZER

Refrigerators

Model	Cu. Ft.	Defrost	Door	Shelves	W"	D"	H"	RL	H.P.	Amps
FR051WWW/0	5	Cycle	1 Solid	3 Adj.	23 3/8	24 1/4	32 3/16	R600a	1/9	1.2

Freezer

Model	Cu. Ft.	Defrost	Door	Shelves	W"	D"	H"	RL	H.P.	Amps
FF051WWW/0M	4	Manual	1 Solid	2 Fixed	23 3/8	24 1/4	32 3/16	R600a	1/8	1.4

GENERAL PURPOSE UNDERCOUNTER REFRIGERATOR AND FREEZER

Refrigerators

Model	Cu. Ft.	Defrost	Door	Shelves	W"	D"	H"	RL	H.P.	Amps
PF051WWW/0	5	Cycle	1 Solid	2 Adj./ 1 Fixed	20 1/4	21 3/4	32 1/4	R600a	1/7	.90

Freezer

Model	Cu. Ft.	Defrost	Door	Shelves	W"	D"	H"	RL	H.P.	Amps
PF041WWW/0M	4	Manual	1 Solid	2 Fixed	21 1/2	22 1/4	32 3/4	R600a	1/7	.90



NEW
LABORATORY UNDERCOUNTER
REFRIGERATORS &
FREEZERS

WHITE DIAMOND SERIES BUILT-IN UNDERCOUNTER REFRIGERATORS AND FREEZERS

The White Diamond Series Built-In Undercounter Refrigerators and Freezers feature a microprocessor temperature controller and display, which allows for precise temperature control, verification and recovery versus mechanical style thermostats. An interior fan driven, forced-draft air circulation system provides excellent temperature uniformity and recovery after door openings (Refrigerators only). Critical samples and supplies are security backed with an extensive package of audible and visual alarms.

****The Undercounter Built-In models are designed with a "front-breathing" ventilation system, eliminating the need for clearance around the top and sides of the unit; enabling the user to incorporate these units into cabinets and other applications where space is limited.**



Environmentally friendly, natural hydrocarbon (HC) refrigerants vastly reduce global-warming potential while meeting new EPA/SNAP mandates.



PR051SS5/0

Features:

- Microprocessor temperature controller
- Digital temperature display, °C/°F switchable
- Audible and visual high and low temperature alarms
- Remote alarm contacts 2 Probes; 1 in air and 1 in sample simulated bottle
- Adjustable operating temperature range: 2°C - 8°C (6°F - 46°F)
- Keyed door lock
- Back wall evaporator cover (select models)
- Leveling legs
- Forced draft air circulation
- Magnetic door gasket for positive seal
- Probe access port (3/8")
- High density urethane foam cabinet and door insulation
- HFC - Free refrigerant (R600a)
- ADA compliant units available
- Three year parts and labor warranty
- UL/C-UL listed
- **In stock with next day shipping**

WHITE DIAMOND SERIES FREESTANDING UNDERCOUNTER REFRIGERATORS AND FREEZERS

The White Diamond Series Freestanding Undercounter Refrigerators and Freezers feature a microprocessor temperature controller and display, which allows for precise temperature control, verification and recovery versus mechanical style thermostats. An interior fan driven, forced-draft air circulation system provides excellent temperature uniformity and recovery after door openings. (Refrigerators Only) Critical samples and supplies are security backed with an extensive package of audible and visual alarms.

****NOTE: Freestanding Models require 4" of clearance around the back and sides of the unit for proper ventilation.**



Environmentally friendly, natural hydrocarbon (HC) refrigerants vastly reduce global-warming potential while meeting new EPA/SNAP mandates.



PR061WWG/0

Features:

- Microprocessor temperature controller
- Adjustable operating temperature range: 2°C - 8°C (6°F - 46°F)
- Digital temperature display, °C/°F switchable
- 2 probes; 1 in air and 1 in sample simulated bottle
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- Magnetic door gasket for positive seal
- Keyed door lock
- Back wall evaporator cover (Select models)
- Forced draft air circulation
- Levelling legs
- Probe access port (3/8")
- High density urethane foam cabinet and door insulation
- HFC - Free refrigerant (R600a)
- Three year parts and labor warranty
- Defrost: Cycle defrost
- UL/C-UL listed
- **In stock with next day shipping**

WHITE DIAMOND SERIES FLAMMABLE STORAGE UNDERCOUNTER REFRIGERATORS AND FREEZERS

In addition to all the features of the White Diamond Series Freestanding Undercounter Refrigerators and Freezers; The White Diamond Series Flammable Storage Undercounter Refrigerator and Freezer meet NFPA and OSHA guidelines 45 and 70 for the refrigerated storage of flammable (volatile) materials in general laboratory work areas. There are no electrical components located in the interior of the unit and the compressor's electrical components have been sealed in a vapor-proof enclosure for additional safety. **These units are NOT designed for use in environments where volatile / explosive conditions could potentially exist.**

****NOTE: Freestanding Models require 4" of clearance around the back and sides of the unit for proper ventilation.**



Environmentally friendly, natural hydrocarbon (HC) refrigerants vastly reduce global-warming potential while meeting new EPA/SNAP mandates.



FF051WWW/0M

Features:

- Microprocessor temperature controller
- Digital temperature display, °C/°F switchable
- Audible and visual high and low temperature alarms
- Remote alarm contacts 2 Probes; 1 in air and 1 in sample simulated bottle
- Adjustable operating temperature range: 2°C - 8°C (6°F - 46°F)
- Keyed door lock
- Back wall evaporator cover (Select models)
- Levelling legs
- Forced draft air circulation
- Magnetic door gasket for positive seal
- Probe access port (3/8")
- High density urethane foam cabinet and door insulation
- HFC - Free refrigerant (R600a)
- ADA compliant units available
- Three year parts and labor warranty
- UL/C-UL listed
- **In stock with next day shipping**

GENERAL PURPOSE UNDERCOUNTER REFRIGERATORS AND FREEZERS

The General Purpose Freestanding Undercounter Refrigerator and Freezer are intended for general purpose applications only where precise temperature control and maintenance are not essential and there are no requirements for temperature alarms. The General Purpose Freestanding Undercounter Freezer features a manual defrost operating system to limit product-damaging temperature spikes found in auto defrost units. Periodic manual defrosting is required.

****NOTE: Freestanding Models require 4" of clearance around the back and sides of the unit for proper ventilation.**



Environmentally friendly, natural hydrocarbon (HC) refrigerants vastly reduce global-warming potential while meeting new EPA/SNAP mandates.



PR041WWW/0M

Features:

- Mechanical thermostat
- Adjustable operating temperature range: 36°F - 46°F (2°C - 8°C)
- Recessed door handle
- Probe access port (3/8")
- White cabinet door
- Cycle defrost refrigerator, manual defrost freezer
- Magnetic door gasket for positive seal
- Levelling legs
- High density urethane foam cabinet and door insulation
- Keyed door lock
- Magnetic door gasket for positive seal
- HFC - Free refrigerant (R600a)
- One year parts and labor warranty
- UL/C-UL listed
- **In stock with next day shipping**