

IL2RB Polyclonal Antibody

Catalog number: 13602-1-AP


Size: 22 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

IL2RB; CD122; P70-75



HEK-293 cells were subjected to SDS PAGE followed by western blot with 13602-1-AP(IL2RB Antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Positive WB detected in:	HEK-293 cells
Calculated IL2RB MW:	551aa,61kd
Observed IL2RB MW:	70-75kd
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3929
GenBank accession number:	BC025691
Gene ID (NCBI):	3560
Full name:	Interleukin 2 receptor, beta

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.