



# Intelli Mixer RM-2S, Small



## Product Description

Powered by V-Spin – the world's most advanced, patented mixing and liquid-handling technology – the Intelli-Mixer RM-2 is truly the latest in cutting-edge science. V-spin harnesses the energy of an unequaled stepper motor, allowing for a wide variety of mixing, vibrating, and rotating cycles. This multipurpose unit is a rotator, vortex mixer, and blotter all in one. Designed to be cold room safe, and with an unprecedented 34 modes of operation, this versatile piece is unrivaled in today's market.

## Features

- Rocker, Vortex Mixer, and Blotter
- 34 automated intelligent mixing programs
- 3 Year Warranty

## Specifications

Rotation Speed	1 - 99 RPM
Rotation Speed Settings Step	1 RPM
Overall Number of Programs	28 Standard + 3 Custom
Environmental Temperature	10 - 45 °C
Relative Air Humidity (at 20 °C)	up to 80 %
Power Supply Adapter	DC 12V, 1.5 A

## Accessories



### IMIX-24

Standard Racks  
Tube Size: 9 mm  
Number of Tubes: 26



### IMIX-23

Standard Racks  
Tube Size: 11 mm  
Number of Tubes: 28



### IMIX-22

Standard Racks  
Tube Size: 13 mm  
Number of Tubes: 24



### IMIX-21

Standard Racks  
Tube Size: 16 mm  
Number of Tubes: 14



### IMIX-20

Standard Racks  
Tube Size: 30 mm  
Number of Tubes: 8



### IMIX-25

Mix Racks  
# of 11mm Tubes: 2  
# of 13mm Tubes: 4  
# of 16mm Tubes: 4  
# of 30mm Tubes: 2



### IMIX-27

Rocking Platform



### IMIX-26

Blotting Rack  
Tube Size: 50 ml  
Holds 6 Tubes



Order: (747) 226-3330

info@elminorthamerica.com

5115 Douglas Fir Rd Suite G | Calabasas, California 91302

elminorthamerica.com

© ELMI North America, Inc. August 2016