



DWK LIFE SCIENCES  
DURAN WHEATON KIMBLE

SOCOREX® QUALITIX®  
Sterile Tips

**BUY ONE  
AND GET  
ONE FREE**

Offer valid until December 31, 2024



Boost your lab's efficiency –  
for every pipette tip you buy,  
get another one for free.

Qualitix® Pipette Tips are all about quality and compatibility. They're free of detectable DNase, RNase, Pyrogen, and ATP, with filter tips certified free of human DNA and protease. Available in rack-style, with or without filters, and low protein binding resin, they offer accurate liquid delivery. Compatible with all major brands, including Socorex®.

Stock up now and double your supply at no extra cost. **Don't miss this limited-time offer!**



customerservice@thomassci.com



833.544.7447



thomassci.com



# SOCOREX® QUALITIX® Sterile Tips – Advantages

## LOW RETENTION – MAXIMAL RECOVERY

Low binding tips are made of silicone-free resin, greatly reducing surface adhesion and providing maximal liquid recovery. Available in 10, 200 and 1000 µL capacities, in sterile and non-sterile racks. Also protease-free certified, they are best suited for proteomic research and accurate handling of viscous liquids.

## AEROSOL PROTECTION – NO CONTAMINATION

Filter tips include a 10 µ pore hydrophobic HDPE filter to retain any contaminant, efficiently protecting both pipette and sample. The high tip purity, as well as the gamma irradiation are certified.



## ORDERING INFORMATION – SOCOREX® QUALITIX® STERILE TIPS

Thomas Number	Cat. Number	Product Description	Case
20A00U420	3090.0010FRS	10µL Filter Fine Point Graduated Tips	12x96 / rack
20A00U423	3090.0010PRS	10µL Low Bind Fine Point Graduated Tips	12x96 / rack
20A00U439	3080.0020FRS	20µL Filter Bevel Point Tips	10x96 / rack
20A00U440	3080.0100FRS	100µL Filter Bevel Point Tips	10x96 / rack
20A00U433	3070.0200FRS	200µL Filter Bevel Point Graduated Tips	10x96 / rack
20A00U436	3070.0200PRS	200µL Low Bind Bevel Point Graduated Tips	10x96 / rack
20A00U428	3070.1000FRS	1000µL Filter Fine Point Tips	6x96 / rack
20A00U438	3070.1000PRS	1000µL Low Bind Fine Point Graduated Tips	6x96 / rack

## To claim your redemption (redeem by 1/31/2025)

1. Obtain your Proof of Purchase (dated between 6/1/24 - 12/31/24).  
Packing slip, purchase order and invoice are all considered valid for proof of purchase.
2. Visit <https://www.dwk.com/na/dwkdeals>
3. Attach your Proof of Purchase and fill-out the online form.
4. The value of freeware will be calculated off of DWK list price and units of measure found on [www.dwk.com](http://www.dwk.com). Free items must be of equal or lesser value. Not for resale.

Offer valid for US and Canadian customers only. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. Offer valid while supplies last. Not valid on prior purchases. Void in government accounts or where prohibited by law or company policy. Prices subject to change without notice. Trademarks contained herein are the property of their respective owners.

001 / EN-06-2024



Thomas  
Scientific

ISO  
9001:2015

833.544.7447

customerservice@thomasci.com

thomasci.com



© 2024 DWK Life Sciences

The trademarks used are owned by DWK Life Sciences GmbH,

DWK Life Sciences LLC or DWK Life Sciences Ltd.

Socorex® and Qualitix® are registered trademarks from Socorex Isba SA.