



Clarity Drug Testing Cups







Screw lid for a no-leak cup.

Easy-to-use.

Generic product available.


510K cleared and CLIA Waived.

Drugs of Abuse	Abbreviation	Cutoff
Amphetamine	AMP	1000 ng/ml
Barbiturates	BAR	300 ng/ml
Buprenorphine/Norbuprenorphine	BUP/NBUP	10 ng/ml
Benzodiazepine	BZO/BZD	300 ng/ml
Cocaine	coc	300 ng/ml
Methamphetamine	mAMP/MET	1000 ng/ml
MDMA or Ecstasy	MDMA	500 ng/ml
Morphine	МОР	300 ng/ml
Methadone	MTD	300 ng/ml
Opiate	OPI	2000 ng/ml
Oxycodone	ОХҮ	100 ng/ml
Phencyclidine	PCP	25 ng/ml
Propoxyphene	PPX	300 ng/ml
Tricyclics	TCA	1000 ng/ml
Marijuana/Hashish	THC	50 ng/ml

Read the drug test results at 5 minutes, do not interpret the results after 10 minutes.







Clear Vision - Clear Results

Clarity Brand Rapid Diagnostics

