

FTSJ1 Polyclonal Antibody

Catalog number: 11620-1-AP


Size: 43 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

FTSJ1; CDLIV, FtsJ homolog 1
(E. coli), FTSJ1, JM23, MRX44,
MRX9, Protein ftsJ homolog 1,
SPB1, TRM7



HepG2 cells were subjected to SDS PAGE followed by western blot with 11620-1-AP (FTSJ1 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated FTSJ1 MW:	329aa, 36 kDa
Observed FTSJ1 MW:	36 kDa
Positive WB detected in:	HepG2 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2199
GenBank accession number:	BC023584
Gene ID (NCBI):	24140
Full name:	FtsJ homolog 1 (E. coli)

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.