

## PAGE2 Polyclonal Antibody

Catalog number: 11961-1-AP


Size: 55 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

PAGE2; G antigen family E member 2, GAGEC2, GAGEE2, PAGE 2, PAGE2



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using 11961-1-AP (PAGE2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate hyperplasia tissue slide using 11961-1-AP (PAGE2 Antibody) at dilution of 1:50 (under 40x lens)

### Applications

Tested applications:	ELISA, IHC
Species specificity:	Human; other species not tested.
Calculated PAGE2 MW:	12 kDa
Observed PAGE2 MW:	
Positive IHC detected in	Human prostate hyperplasia tissue
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag2566
GenBank accession number:	BC054022
Gene ID (NCBI):	203569
Full name:	P antigen family, member 2 (prostate associated)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.